Pink tens

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5 passage Saint Bernard 75011 Paris, France.

Bluetens limited 

Unit G, 15/F, TAL Building, 49 Austin Road, Kowloon, Hong Kong
Contact : cs@bluetens.com

S4M EUROPE SAS 

59, rue Castellion 01100 Oyonnax FRANCE - regulation@s4m-europe.com

English	4
Français	18
Deutsch	32
Español	46
Other languages	61
Electrode placements	64
Description of currents output	62
Electromagnetic compatibility	66

Welcome!

Pinktens is a medically certified TENS device, designed to support women with all perineal-related concerns. Please read this user manual carefully for safe and optimal use of the device.

WHO ARE WE?

For the past 10 years, we have been developing electrostimulation devices to address various medical conditions, particularly those related to menstrual and endometriosis pain. Pinktens was developed by pelvic floor specialists and tested in hospitals on a wide range of patients suffering from urinary or fecal incontinence.

We were committed to helping all affected women regain well-being and confidence in their daily lives. From now on, Pinktens is taking care of your perineal health :)

This user manual is available by following the link:
<https://www.pinktens.com/en/content/48-user-manual>

1


2

3


4

5

6

- 1 Pinktens (pink) CL2P (x1)
- 2 White USB-A to USB-C cable, 80 cm (x1, PN: V2C-A1C080WH)
- 3 Electrode cable, single output, 25 cm (x1, PN: V2E-C1S025WH)
Electrode cable, dual output, 110 cm (x1, PN: V2E-C2S110WH)
- 4 Mixed electrode pack: 4S (50x50 mm) + 4M (100x50 mm), PN: MG-MIXP
- 5 User manual (x1)
- 6 Protective pouch (x1)

EN PRECAUTIONS FOR USE

Target patient population: (Adults: 18 years old and above)

TENS Function Indications: Patients suffering from acute or chronic pain (e.g., neck, shoulder, back, or leg pain), pelvic pain, primary/secondary dysmenorrhea, or endometriosis-related pain (with confirmed pain etiology and medical clearance to ensure absence of contraindications).

Eligible Operator Population:

- Must be 18 years old or above, with full civil capacity;
- If the patient operates the device themselves: Must have basic reading and comprehension skills to understand the instruction manual and follow steps, and be able to independently report physical sensations (e.g., pain, numbness) during use;
- If operated by medical staff/family members: Must be familiar with device functions and operating procedures, have basic safety knowledge of electrical stimulation, be able to assist patients in adjusting device parameters, and monitor the patient's physical state during use.

Contraindicated Operator Population

- Minors under 18 years old;
- Individuals with visual, auditory, or physical impairments that prevent accurate device operation;
- Individuals lacking basic safety awareness and unable to identify abnormal conditions during patient use

MAIN DEVICE:

- The electrode is treated as a BF-type applied part.
- Device name: Pinktens
- Product model: CL2P
- Weight: 33 grams
- Battery voltage: 3.7V 500mAh lithium-ion polymer battery
- USB charging: Type C 5V/500mA
- Frequency: 1Hz to 1200Hz
- Output pulse: 24 to 400uS
- Output voltage: 0V to 120V (at 500Ω load)
- Output power: 10mW (r.m.s) MAX
- Operating environment:
 - Temperature: 5°C to 40°C
 - Humidity: 30%RH to 80%RH
 - Atmospheric pressure: 700 to 1060hPa
- Storage environment:
 - Storage temperature: -20°C to 70°C
 - Humidity: 30%RH to 80%RH
 - Atmospheric pressure: 700 to 1060hPa
- Product dimensions: 50.8mm X 51.5mm X 18mm
- Minimum lifespan: 5 years
- Software version: V3.3

ELECTRODES:

- Support material: non-woven fabric
- Safety: Biocompatibility of the G607 gel used in this product is compliant with ISO 10993-5, ISO 10993-10 standards
- Lifespan: sealed storage for 3 years, to be stored at room or cool temperature (5°C to 35°C) and protected from direct sunlight.

CONTRAINDICATIONS

1. Do not use this device if you have a pacemaker, implanted defibrillator, or any other implanted metallic or electronic device. Such use could cause electric shock, burns, electrical interference, or death.
2. The device should not be used when cancerous or other lesions are present in the treatment area.
3. Stimulation should not be applied over open wounds, rashes, or swollen, red, infected, or inflamed areas (e.g., phlebitis, thrombophlebitis, varicose veins, obliterative arteriosclerosis, etc.).
4. Electrode placements should be avoided in the carotid sinus region (front of the neck) or transcranially (across the head).
5. Anxious users - the use of electrical stimulation requires user cooperation; therefore, the procedure should not be attempted in users with communication disabilities or mental incapacity.
6. Users with cerebrovascular issues - individuals with a history of aneurysm, stroke, and transient ischemia should not undergo electrical stimulation, as it stimulates peripheral blood flow, which could be fatal in such cases.
7. Epileptic users - Electrical stimulation impulses have the potential to trigger a seizure.
8. Cases of acute pain/pain of unknown origin - using an electrical stimulator in undiagnosed cases may hinder diagnosis.
9. Acute inflammation and swelling of the genitourinary system, tumors, urinary stones, history of nervous system diseases, and other illnesses.
10. Do not use during pregnancy, especially during the first trimester.
11. Presence of a severe systemic disease, such as scleroderma, diabetes mellitus with severe complications, myopathy, multiple sclerosis, Parkinson's disease, myelopathy, or serious diseases of the liver, lungs, kidneys, hematological system, or other comorbidities.

WARNINGS

The use of the device should be discontinued in the following circumstances:

1. If you feel unwell or uncomfortable, stop immediately. If there are any abnormalities in the skin or body, consult your doctor and follow their instructions.
2. If you wish to move the electrode to other parts of the body during operation, please turn off the power before moving it. Failure to do so may result in an electric shock.
3. Do not use this device while wearing other electronic devices.
This could cause malfunction of these devices (e.g., malfunction of an electronic watch).

4. As subcutaneous use presents a high risk to humans or other devices, the use of this device is prohibited.
5. Use under supervision mandatory for children (over 7 years old) and disabled persons. Risk of accident or discomfort.
6. Do not use this appliance in a bathroom or a place with high humidity. Do not use while bathing or showering. Strong electric shock could occur.
7. Do not use while sleeping. This could damage the product, and the displacement of the electrode could cause discomfort.
8. Do not use while driving. It could cause an accident due to the strong stimulation.
9. Do not keep the electrodes on your body all the time. This could cause skin irritation.
10. Do not use metal to touch the electrode. This could cause an electric shock.
11. Application near the thoracic region could cause cardiac discomfort. Seek the advice of a doctor beforehand.
12. If in doubt, consult your doctor.
13. Discontinue use and do not increase the intensity level if you experience discomfort during use.
14. The electrode cable must not encircle your neck to avoid strangulation.
15. Do not modify the equipment. Please use original accessories, it will cause an unusable effect if you use third party accessories on the market.
16. Keep out of reach of children to avoid inhalation or ingestion of small parts.

PRECAUTIONS FOR USE

1. If you have a skin allergy, you should stop using it immediately.
2. Stimulation must not be applied through or over the head, directly to the eyes, covering the mouth or on the front of the neck (especially the carotid sinus).


SIDE EFFECTS

1. Users may experience skin irritation and a burning sensation under the stimulation electrodes applied to the skin.
2. Users may experience pain or discomfort during or after application.

Note: Users should stop using the device and consult their doctor if they experience any adverse reactions.

THE DEVICE

The Pinktens turns on and off by pressing the blue button in the center of the device. Here is an explanation of the inputs, outputs, and buttons of the device.

CHARGING THE BATTERY

The Pinktens is charged using the USB-C cable included in the box.

Charging time is approximately 1 hour.

When the battery is below 20%, the charge indicator light flashes orange.

This indicator light will remain orange during the charging time and will go out once the battery is full.

We recommend charging the device with its USB-C cable connected to a computer.

If you wish to use a mains adapter, please ensure that the power is compatible.

The battery recharges on average once a month, if the device is used once every 2 days.


Note:

- This device cannot be used while charging.
- The mains adapter must comply with standard IEC62368-1, respect at least 2 MOOP and have an output of 5V/1A.


CONNECTING THE ELECTRODES

Here's how to connect the electrodes:

1. Short electrode cable + 2 electrodes (S or M)

2. Long electrode cable + 2 electrodes (S or M)

DOWNLOAD AND INSTALLATION

The Pinktens Application is available on Google Play and Appstore.

The Application is compatible with Smartphones and Tablets:

iPhone : iOS 13 and above

Android : Android 10 and above

(Modified Android versions may make the Application incompatible)


Download and install the Application according to your Smartphone.

Notes:

Bluetooth must be activated to use the Application.

Android smartphones must also authorise access to the location service, to the Application.

The Pinktens Application does not have access to any data. This condition is required by Android versions to stabilise Bluetooth.

OVERVIEW AND USAGE

> Choose from the main functions on the home screen:

ELECTROSTIMULATION PROGRAMS

15 programs for pain relief

TREATMENT PROTOCOLS

Incontinence, postpartum, and prolapse therapies

WELLNESS PROTOCOLS

Prevention and sexual comfort

MAIN DEVICE


- Remove the cables from the device after use.
- Use a soft cloth to clean dirt on the device
- Use a damp cloth with a neutral detergent or alcohol to clean if necessary.
- Do not use thinner or petrol to clean.
- Do not use alcohol on the back of the device
- Do not use any other detergents.
- Keep water or liquids away from the device
- Keep out of direct sunlight, heat, damp and dusty places.

ELECTRODES

- Do not use extra gel with the electrodes.
- Always replace and store the electrodes on the correct side of the plastic provided.
- Use a little water to clean the electrode if necessary.
- Do not wash.
- Do not wipe the electrode with a cloth or paper.
- The lifespan of the electrodes varies between 30 and 50 uses depending on the quality of the skin, hairiness, perspiration and the care taken.
- Electrodes are consumables; when they no longer adhere to the skin, it's a sign that they need changing.
- Do not stack or bend the electrodes.
- Do not touch the gel with your fingers.
- Apply to clean, dry skin, without cream or oil.

The product must not be disposed of with household waste.
Consumers are legally obliged to dispose of electrical equipment at a collection point in their town.
Contact your local council for more information or a retailer.

Essential informations for a correct use is indicated by the corresponding symbols. The following symbols can be seen on the device and on the labelling.

SYMBOL	TITLE
IP22	IP code of the device: the degree of protection of this device against the intrusion of solid foreign objects -- ≤ 12.5 mm in diameter (and against access to dangerous parts with a finger); the degree of impermeability is of the drip type (inclination of 150).
	Disposal in compliance with directive 2002/96/EC (WEEE)
	This side up
	Date of manufacture
	Batch code
	CE mark, Certificate issued by TUV SUD.
	Manufacturer
	Type BF applied part
	Warning: to avoid the risk of injury, read the manufacturer's instructions carefully.
	Symbol for «AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY»
	Fragile, handle with care
	Store the product in a dry place. Protect from water and rain.
	Product packaging is recyclable
	Importer
	Medical device

DO YOU NEED HELP?

Our customer service team deals directly with all requests for advice or after-sales service. Our agents are fast and available by e-mail at support@pinktens.com.

We'll get back to you within a day (excluding weekends and public holidays).

Online troubleshooting: <https://www.pinktens.com/fr/content/28-depannage>

The Pinktens is guaranteed for two years. This guarantee does not cover damage caused by improper maintenance or use that does not comply with the relevant documentation, accidents (impact, fire, liquid), modifications, repairs, misuse or normal wear and tear. This warranty covers defects in materials or workmanship. This warranty only covers products purchased directly from Bluetens or an authorized Bluetens reseller.

In the event of a defect covered by the warranty, Bluetens will send you a return slip and will be responsible for the repair or exchange of your device as well as its return.

WARRANTY PERIOD

Your device is guaranteed for 2 years from the date of receipt of the product if purchased online, or from the date of purchase if purchased from a retailer.

NOT COVERED BY THE WARRANTY

Electrodes, cables and accessories included in the box are not covered by this warranty.

WARRANTY DISCLAIMER

Bluetens has not authorised any agency to repair its products. If you experience any problems with your device, please contact us. Bluetens will not be responsible for the results of any repair. The intervention of a third party on our products will void the 2-year warranty. The use of third-party electrodes or accessories could cause the device to malfunction and will void the warranty. The user should not attempt to repair the device or its accessories by himself. Each device has been inspected through systematic validation. The performance is stable and does not require calibration or validation. If your product does not achieve the expected performance and the basic function changes during normal use, please contact us. Please note that Bluetens will not be responsible for quality issues caused by improper use.

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
CHEVILLE
KNÖCHEL
TOBILLO

DESCRIPTION OF CURRENTS OUTPUT

Instant electric voltage/current when using Electrodes:

Load 500Ω: MAX 120V / MAX 240mA

Load 1000Ω: MAX 150V / MAX 150mA

TENS

N°	PROGRAM	OUTPUT PARAMETERS
5	Endorphin 1	Frequency: 5Hz - Impulse: 250μs - Treatment time: 30 minutes
6	Chronic pain 1	Frequency: 80Hz - Impulse: 200μs - Treatment time: 60~120 minutes

OTHER LANGUAGES

**THIS USER MANUAL IS AVAILABLE IN OTHER LANGUAGES
AT THE FOLLOWING LINK:**

um.pinktens.com

ELECTROMAGNETIC COMPATIBILITY

The device has been thoroughly tested and inspected to assure proper performance and operation. The EMC Test Report Number of TUV: 68.730.20.0074.01. Model: CL2P

TEST CONFIGURATION

Configuration	Description
AC Powered	230VAC/50Hz
AC Powered	100VAC/50Hz & 240VAC/50Hz (Only for Voltage dips and Voltage interruptions)

OUTPUT PARAMETERS

Output Frequency	1Hz	10Hz	50Hz	100Hz	200Hz	500Hz	1KHz	1,2KHz
Output Pulse	400uS	300uS	100uS	100uS	50uS	50uS	34uS	24uS
Output Voltage V _{pp}	0-120V	0-120V	0-110V	0-85V	0-85V	0-57V	0-57V	0-57V

MODES OF OPERATION:

Mode	Description
Charging	The EUT's USB-C port is connected to a single-port USB-C cable (length: 80 cm) and is charged by an AC / DC adapter.
TENS + Lien BT	The EUT is powered by a rechargeable 3.7VDC lithium battery. Its USB-C port is connected to a dual-electrode cable (4 ends, length: 110 cm), and the electrode cable terminals are connected to electrode pads (Size: 50 x 100 mm). The part of the electrode cable that is 40 cm away from the EUT is placed in 0.9% saline solution. The EUT is connected to the APP via Bluetooth, and its processing functions are controlled by the APP or by itself. Next, the EUT is set to TENS mode and adjusted to maximum power for operation.

EMISSION TESTS EN 60601-1-2:2015, IEC 60601-1-2:2014

Test Condition: Group1, Class B	Test Result
Radiated Emission 30MHz to 1000MHz	Pass
Conducted Emission on AC 150kHz to 30MHz	Pass
Harmonic Class A	Pass
Flicker	Pass

IMMUNITY TESTS EN 60601-1-2:2015, IEC 60601-1-2:2014

Test Condition: Group1, Class B	Test Result
Electrostatic Discharge (IEC 61000-4-2) ±2kV ±4kV ±6kV ±8kV ±15kV	Contact: ±8kV Air: ±2kV ±4kV ±8kV ±15kV
Radiated Immunity (IEC 61000-4-3) 80MHz to 2700MHz 10V/m 385MHz 27V/m 450MHz 28V/m 710MHz, 745MHz, 780MHz 9V/m 810MHz, 870MHz, 930MHz 28V/m 1720MHz, 1845MHz, 1970MHz 28V/m 2450MHz 28V/m 5240MHz, 5500MHz, 5785MHz 9V/m	10V/m, 80% Am at 1kHz 27V/m PM at 18Hz 28V/m FM ± 5 kHz deviation at 1kHz sine 9V/m PM at 217 Hz 28V/m PM at 18Hz 28V/m PM at 217 Hz 28V/m PM at 217 Hz 9V/m PM at 217 Hz
Electrical Fast Transient (IEC 61000-4-4) For AC Power port: ±2kV, 100 kHz repetition frequency For signal input/output port: ±1kV, 100 kHz repetition frequency	For AC port: Power supply lines: ±2kV
Surge (IEC 61000-4-5) For AC Power port: ±0.5kV, ±1kV (Line to line) ±0.5kV, ±1kV, ±2kV (Line to Ground) For signal input/output port: ±2kV (Line to Ground)	For AC port: (L-N): ±0.5kV ±1kV
Conducted Immunity (IEC 61000-4-6) For AC Power port, patient coupling port and signal input/output port: 150KHz to 80MHz 3V (80% Am at 1kHz) ISM and amateur radio bands between 150KHz to 80MHz 6V (80% Am at 1kHz)	For AC port and patient coupling port: 150KHz to 80MHz 3V (80% Am at 1kHz) ISM and amateur radio bands between 150KHz to 80MHz 6V (80% Am at 1kHz)
Voltage Dips and Interruption (IEC 61000-4-11) 0%, 70%, 0% of UT	For AC power port: 0% for 0.5 cycle (at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°), 0% for 1 cycle (at 0°), 70% for 25 cycles (at 0°), 0% for 250 cycles
Power Frequency Magnetic Field (IEC 61000-4-8) 50Hz, 60Hz 30A/m	50Hz: 30A/m 60Hz: 30A/m

Take care of yourself

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Cuidese

pinktens.com